

## Rigaku MiniFlex Standard Operating Procedure

Manager: Dr. Dmitri Barbash  
Senior Research Scientist


Contact: Drexel University  
Bossone 112  
3141 Chestnut Street  
Philadelphia, PA 19104  
Office: 215-571-4611  
Mobile: 267-918-3188  
[dvb36@drexel.edu](mailto:dvb36@drexel.edu)

Latest Update: 6/24/2020

*These notes are meant to serve as an aid to assist users who have been trained and certified by MCC Staff. If ever you are unsure about the correct operation of the instrument or any of its components, please consult a MCC staff member before continuing. Never use equipment that you are not trained and approved to use.*

*Before using the MCC, please review the MCC User Handbook available through our website.*

Step 1. Open up the program that operates MiniFlex. It is called MiniFlex Guidance.  
Skip to the Step #3 if the MiniFlex is already running




Step 2:

Log in Name : administrator

Password is blank, aka no password required




Step 3. MiniFlex operating program opens up in the following format.



There are 2 useful commands that can be performed anytime.

First. You can Initialize Layout under the View. This clears out previous user window arrangements.

Second. Make sure you are in the Package Measurement under Tasks Menu.



Step 4. Click the yellow Package Measurement to proceed with the scan set up.



General measurement Window is opened up at this point.


Automatic generation of data file names

No.	Exec.	Folder	File name	Sample name	Memo	Condition	Analysis condition	Spin	Sample
1	<input checked="" type="checkbox"/>	C:\XRD data\Arkita\02072020\'	DMAC.ras	...		No1	Set...	<input type="checkbox"/>	1
2	<input type="checkbox"/>			...		No1	Set...	<input type="checkbox"/>	1
3	<input type="checkbox"/>			...		No1	Set...	<input type="checkbox"/>	1
4	<input type="checkbox"/>			...		No1	Set...	<input type="checkbox"/>	1
5	<input type="checkbox"/>			...		No1	Set...	<input type="checkbox"/>	1
6	<input type="checkbox"/>			...		No1	Set...	<input type="checkbox"/>	1
7	<input type="checkbox"/>			...		No1	Set...	<input type="checkbox"/>	1
8	<input type="checkbox"/>			...		No1	Set...	<input type="checkbox"/>	1
9	<input type="checkbox"/>			...		No1	Set...	<input type="checkbox"/>	1
10	<input type="checkbox"/>			...		No1	Set...	<input type="checkbox"/>	1

XG termination condition:     Optical devices confirmation message:

The status quo     Stop the X-Ray     Show     Don't show

Press the Measurement conditions to set up the parameters for your scans. All physical meaning of the parameters are discussed during training.



Step 5 (optional). X-Ray turning on Procedure is under XG control and 40Kv, 15mAmp is the max values for the Xray power. Xray usually stays on 24/7.

Locate the feedback window and make sure Tube voltage and tube Current is reading zero. Proceed to XG control, turn on Generator and Bring up Tube voltage and Current to 40kV and 15mAmp.

